EN

ENAMELITE SPRAY GLAZE INSTRUCTIONS

Indications: Enamelite Ceramic Spray Glazes are intended to be used for aerosol applications of ceramic porcelain glazes for dental restorations (crowns and bridges).

- Shake well until mixing ball moves freely inside can.
- Prepare restoration(s) in usual manner for stain and glaze.
- Use glaze on low fusing PFM porcelain, pressed ceramic, Lithium Disilicate and Zirconia restoration(s).
- Hold spray can upright 2-4 inches away from restoration(s) over dried and unfired stains.
- Use short bursts while shaking can between sprays on buccal, lingual, mesial, distal and occlusal surfaces
- If renewed application is necessary, spray restoration in the same manner again.
- If multiple units are placed on ceramic pillow, spray units simultaneously while shaking can often.

Firing Schedule:

Pre-Dry: 3-5 minutes High Temp: 1346°F(730°C) -1562°F(850°C)
 Start Temperature: 572°F/300°C Vacuum: None
 Firing Rate: 72°F/40°C Cool Time in Firing Schedule: None
 CTE: 10.0 ± .5 X 10⁻⁶/°C Hold Time: 1 minute

NOTE: If glaze should settle on the underside of restoration(s), dust with small bristle brush. Clear valve after each application by turning can upside down and depressing spray button.

WARNING: Do not inhale aerosolized ceramic dust. Do not get dust or liquid in contact with eyes or skin. Personal Protective equipment such as gloves, goggles, mask and lab coats are recommended when using this product.

Type I, Class 1b

Refer to product SDS for additional health and safety information.

Keystone Industries
 52 West King St.
 Myerstown, PA 17067 USA
 800.333.3131
 keystoneindustries.com

Keystone Europe LLC
 Boltavenweg 7
 5349BC Oss Netherlands

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