Bosworth New Truliner®

Denture Corrective Relining Material

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Type 2, Class 1 I Autopolymerizable

Powder & Liquid System

Keystone Industries
52 West King St.
Myerstown, PA 17067 USA
800.333.3131
keystoneindustries.com

keystoneindustries.com/symbolglossary

ORDERING INFORMATION
BESTELLINFORMATIONEN
INFORMATIONS CONCERNANT LES COMMANDES
INFORMAZIONI PER LE ORDINAZIONI
INFORMAÇÃO PARA ENCOMENDA
INFORMACIÓN DE PEDIDOS

8 oz (226 g)
Powder

6 oz (177 mL)
Liquid

1/3 oz (10 mL)
Bonding Liquid

1 lb (454 g)
Powder

8 oz (236 mL)
Liquid

1/3 oz (10 mL)
Bonding Liquid


Powder available in pink, vein, and clear.

[REF] 0921972	Standard Kit, Clear
[REF] 0921971	Standard Kit, Vein
[REF] 0921970	Standard Kit, Pink
[REF] 0921979	Bulk Kit, Clear - 1 lb. (454g)
[REF] 0921978	Bulk Kit, Vein - 1 lb. (454g)
[REF] 0921977	Bulk Kit, Pink - 1 lb. (454g)
[REF] 0921984	Bonding Liquid - 4 oz. (118ml)
[REF] 0921986	Liquid Only - 8 oz. (236ml)
[REF] 0921982	Powder Only, Clear - 1 lb. (454g)
[REF] 0921981	Powder Only, Vein - 1 lb. (454g)
[REF] 0921980	Powder Only, Pink - 1 lb. (454g)
[REF] 0921976	Liquid Only - 6 oz. (177ml)
[REF] 0921975	Powder Only, Clear - 8 oz. (226g)
[REF] 0921974	Powder Only, Vein - 8 oz. (226g)
[REF] 0921973	Powder Only, Pink - 8 oz. (226g)


Bosworth New Truliner Instructions for Use

New Truliner is a PEMA denture corrective relining material used to restore a denture's original fit.

DIRECTIONS


1. Clean and dry denture. Roughen denture surface (2mm) and periphery for best adhesion.

2. Coat buccal and labial surfaces where bonding is NOT desired with Bosworth Masque.


3. Using a brush, apply bonding liquid to all surfaces to be relined.


4. Slowly add 1 full plastic vial measure of powder (15ml) to 1 full glass vial measure of liquid (7.6ml) in mixing cup. Spatulate for 30 seconds until material is a honey-like consistency.

5. Apply mixture to all surfaces to be bonded. When material is flat in color and has absorbed all liquid, rinse patient's mouth and insert denture.

6. Have patient close to align occlusion.

7. **Keep denture in patient's mouth until material is in a dough-like stage (approx. 1 minute).**

8. Remove from patient's mouth before final set. Rinse mouth with cold water. Adjust as needed and polish.

NOTES

1. Before making a permanent denture reline, check for sores or cancer and treat the tissue with a tissue conditioner such as Bosworth Softone to eliminate the possibility of the final reline duplicating the disturbed or deformed mucosa caused by the deficient denture.
2. **This product generates heat during polymerization that could be harmful to the patient or dental professional. Immediately remove from mouth when the patient senses the rising heat. Do not cure material to completion in mouth.**
3. Recommended: Bosworth's Pressure Pot for a smoother, color-stable and reduced porous reline.

CAUTION: Certain individuals may be sensitive to products containing acrylic.

MIX RATIO

7.6ml (1 full glass vial) liquid to
15ml (1 full plastic vial) powder

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