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Bosworth Superpaste Instructions for Use

Bosworth Superpaste is a zinc oxide eugenol paste used for preliminary compound impressions, mucostatic impressions, corrective washes, rebasing dentures, and denture relines.

DIRECTIONS

- Extrude equal volumes of Superpaste catalyst (green) and base (white) on mixing pad.¹
- Mix material quickly to uniform color. Approximately 30 seconds to 1 minute.
- Apply mix to tray evenly. Load time is approximately 30 to 60 seconds.
- Insert tray and position slowly without pressure.
- Allow paste to flow to all undercuts and irregularities of the tissue. Sets in less than 3 minutes. Remove from mouth.³
- Rinse with cold water to remove saliva. Disinfect if necessary.⁵

NOTES

- Mix ratio is approximately 13 g base to 5.6 g catalyst. The addition of more base causes the mix to be thicker and take longer to set. More catalyst produces a thinner mix which sets faster.
- Always return the right caps to their respective tubes as an interchange may cause material to set. Keep tubes tightly closed.
- Superpaste is affected by temperature and humidity. Warmer temperatures and higher humidity accelerate set and colder temperatures and lower humidity delay set. Use only at room temperature.
- For a final mucostatic impression in a base plate tray or denture, a corrective wash with Plastopaste can be used.
- To disinfect: Immerse in a 1:213 iodophor for 10 minutes or 2% glutaraldehyde solution for 10 hours, prepared according to the manufacturer's instructions. DO NOT use chlorine compounds.

Type 1: Hard Eugenol (in catalyst): 21.3%
Shelf Life: 2 years