# Bosworth Light Liner<sup>®</sup>

## Light Cure Chairside Reline Composite

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30 cc Air Barrier

#### LIQUID Danger

Causes skin irritation. May cause an allergic skin reaction. Causes serious eye irritation. Harmful if inhaled. May cause respiratory irritation. May damage fertility or the unborn child. Suspected of causing cancer. (Applies to Soft Liquid Only)

#### POWDER Danger

Causes skin irritation. May cause an allergic skin reaction. Causes serious eye irritation. Harmful if inhaled. May cause respiratory irritation. May cause cancer. Suspected of damaging fertility or the unborn child.

#### AIR BARRIER Warning

Harmful if swallowed. Causes eye irritation.

#### Refer to product IFU or SDS for additional health, safety and disposal information.

 This product can expose you to 1,4-Dioxane, which is known to the State of California to cause cancer. This product can expose you to chemicals including cadmium, which is known to the State of California to cause cancer and birth defects or other reproductive harm. This product can expose you to chemicals including Titanium dioxide, which is known to the State of California to cause cancer.


 This product can expose you to chemicals including N,N-dimethyl-p-toluidine, which is known to the State of California to cause cancer, and Di-n-butyl phthalate, which is known to the State of California to cause birth defects or other reproductive harm. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

Keystone Industries  
52 West King St.  
Myerstown, PA 17067 USA  
800.333.3131  
[keystoneindustries.com](http://keystoneindustries.com)

Liquid Flashpoint 230°F (110°C)

[keystoneindustries.com/symbolglossary](http://keystoneindustries.com/symbolglossary)


# Bosworth LightLiner Instructions for Use

Bosworth LightLiner is a unique light-cure chairside reline material that is available in a soft or hard non-methyl methacrylate formula.


## PREPARATION

Remove denture from patient's mouth and grind 1-2mm from the entire palatal surface with a bur to expose fresh acrylic base. Coat patient's teeth and labial surface with a silicone lubricant such as Bosworth Masque to prevent the composite from sticking.


## APPLICATION


1. Add 1 vial of liquid (4g) to plastic cup of powder (6g). Spatulate for 30 seconds to obtain uniform mix. Mix slowly to reduce bubbles.

2. Coat denture with a thin layer of LightLiner mix to serve as a bonding agent.<sup>1</sup>

3. Set aside remainder of mix until it reaches a honey-like consistency (2-3 minutes).

4. Spatulate mix into denture and place in patient's mouth. Wipe off excess. Hold in patient's mouth until composite reaches a rubber-like consistency (3-5 minutes).

## CURING

### Soft Formula

1. Remove denture from patient's mouth, rinse and dry.
2. Place in a light cure unit with a 100-watt halogen for 5 minutes.<sup>2</sup>
3. Remove denture and coat with Air Barrier.
4. Cure for an additional 5 minutes.
5. Wash denture to remove Air Barrier coating.
6. Trim and polish as necessary.

### Hard Formula

1. Remove denture from patient's mouth, rinse and dry.
2. Coat denture with Air Barrier.
3. Place in a light cure unit with a 100-watt halogen for 5 minutes.<sup>2</sup>
4. Wash denture to remove Air Barrier coating.
5. Trim and polish as necessary.

## NOTES FOR USING LIGHTLINER

1. If air bubbles surface, puncture immediately with a spatula to eliminate any porosity.
2. Times listed above are applicable only if curing with a 100 watt oven light. If curing with a 250 watt oven light, time will be approximately 1 1/2 minutes. May be cured with a hand-held unit with a wavelength range greater than 450 nm (approximately 10-15 minutes). A more powerful light source will reduce the time of cure and a less powerful light source will increase the time of cure.

CAUTION: Certain individuals may be sensitive to products containing acrylic.

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