Bosworth Fastray™ LC

Light Cure Custom Tray and Acrylic Base Plate Material

Porte-empreinte individuel et plaque-base en acrylique / Lichtgehärtetes Autopolymerisat für individuelle / Abformlöffel und Kunststoff-Basisplatten / Materiale fotopolimerizzabile per portaimpronte personalizzate e basi provvisorie per protesi acriliche. / Material para placas de base acrílicas e moldadeiras individuais fotopolimerizável / Material Fotocurable para Cubetas Individuales y Bases Acrílicas

Warning

Causes skin irritation. May cause an allergic skin reaction. Causes serious eye irritation. May cause respiratory irritation. **Refer to SDS for additional Health, Safety and Disposal Information.**

CAUTION: US Federal law restricts this device to sale by or on the order of a dental professional.

Flashpoint >150°C (>302°F) / Point d'inflammation >150°C (>302°F) / Flammpunkt >150°C (>302°F) / Punto di infiammabilità >150°C (>302°F) / Ponto de centelha >150°C (>302°F) / Punto de inflamación >150°C (>302°F)

Keystone Industries
52 West King St.
Myerstown, PA 17067 USA
800.333.3131
keystoneindustries.com

keystoneindustries.com/symbolglossary

Made in EU


Fabriqué dans l'U / Herstellt in den EU
Fabbricato in EU / Fabricado nos EU / Fabricado en EU

PR271 BOSLAB-137 Rev 1 08/2021


Bosworth Fastray LC Instructions for Use

Fastray LC is a one-component, pre-formed, light-cured custom tray and acrylic base plate material.


DIRECTIONS


1. Coat prepared model with a thin layer of separating medium.


2. Apply Fastray LC wafer using light pressure with fingers on model.

3. Trim wafer to the size of the model with a pointed instrument.

4. Use excess material to form a handle.


5. Use a blunt instrument to make perforations on the palate.

6. Place in light-cure unit for 1.5-5 minutes depending on light source power.*

7. Remove tray from the model. Clean away any separating material.

8. Turn tray over. Light-cure for 2 minutes. Trim and polish as desired with acrylic bur. Remove inhibitor layer with alcohol wipe.

NOTES

- *1. If curing with a 100 watt oven light, time will be approximately 5 minutes. If curing with a 250 watt oven light, time will be approximately 1 ½ minutes. May be cured with a hand-held unit with a wavelength range greater than 450 nm (approximately 10-15 minutes). A more powerful light source will reduce the time of cure.
2. Fastray LC is water and light sensitive. Keep the box closed. Material may get hard after 15 minutes of light exposure.

CAUTION: Certain individuals may be sensitive to products containing acrylic.

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