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Bosworth Trim/Trim II Instructions for Use

Trim and Trim II are fast-setting, non-irritating PEMA temporary crown and bridge materials.

DIRECTIONS

1. Take an impression with material of choice. (Recommended: Bosworth Supergel)
2. Complete tooth preparations. Take another impression from which the final crown or bridge is to be made. Match tooth shade to closest Trim or Trim II shade.
3. Cut peripheral border and interproximal areas for best fit.
4. Mix powder and liquid until honey-like consistency. Approximate powder to liquid ratio is 6.3g of powder to 3ml of liquid. This amount would make 9 full jacket crowns.
5. Fill each tooth area within the impression where preparatory work has been done with Trim acrylic.
6. When Trim gets dull, reseat the impression to return to its original position.
7. Leave in place for 30-90 seconds. Remove when first feels warm. When material is rubbery it may be removed from impression and trimmed. May be placed back into impression until fully set.
8. Trim and polish as desired. Cement bridge or crown in place with a non-eugenol cement.

NOTES

1. Using Bosworth Gorilla Seal or Glaze saves polishing time. Gorilla Seal will bond directly to the temporary when cured with a halogen bulb.
2. Do not allow material to cure in mouth.
3. May be bench cured or placed in pressure pot (Recommend: Bosworth Pressure Pot).

CAUTION: Certain individuals may be sensitive to products containing acrylic.

WORKING TIME: 3-4 minutes / 2-3 min (fast)
TOTAL SET TIME: 8 minutes / 6 minutes (fast)