EN CLEARMET™

“The Clear Choice”

Indications for Use: Fabrication of removable dental prosthetic devices, such partial dentures, occlusal splints and night guards.

Storage Conditions: Store in a cool dry location in the original packaging until use.

Warnings: Provide adequate ventilation during thermal processing. Do not breathe the vapors from heated material. Respiratory protection, thermal gloves and eye protection are recommended during processing. Allergic reactions to the ClearMet material have not been reported. If the end user experiences discomfort when wearing the finished appliance, they should be instructed to immediately remove it and seek medical attention as needed.

Caution: Federal law restricts this device to sale by or on the order of a dental professional.

Technical Data: Recommended PSI: 105
Recommended Melting Temp: 542°F. Melt time of 16 minutes for small cartridges and 19 minutes for medium cartridges, and a hold time of 2 minutes. **NOTE:** Temperatures may vary depending on machine type.

INSTRUCTIONS FOR USE

ClearMet Framework Injection Procedure

ClearMet is a clear thermoplastic resin specifically designed for durability and biocompatibility. The following procedure outlines how to use ClearMet material with the FlexPress™Digital Injection System. The material is intended for use in an injection molding process.

1. Block out the undercuts on model for correct path of insertion. Place very thin layer of wax over saddle areas on the model. This is in preparation for the acrylic that will be processed later. The thickness used is 1 to 1.5mm thick. Cut a tissue stop into the wax on free end distal extensions. Stability is the reason for the tissue stop.

2. Duplicate Model.

3. Wax up the framework. Be sure when waxing up to evenly distribute the wax. This will help prevent any blockage when injecting the ClearMet material.

4. Coat the inside of the flask with petroleum jelly. This will prevent the stone from adhering to the flask.

5. Invest the framework using regular Buff stone. When the stone sets, place sprues using 3/16th (4.75mm) utility wax. Extend the sprue to the injection port in the flask.

6. Coat the stone with stone separator or petroleum jelly and bolt the top half of the flask in place. The utility wax should fill the injection port of the flask. This will keep the stone from leaking out during the second half of the investing process.

7. Fill up the top section of the flask using Buff Stone pouring and vibrating stone into flask. The stone should not protrude from the widow in the flask.

8. When the top flask section of the stone is set, remove the bolts but leave the flask together and place in boiling water for 4-5 minutes.

9. After all of the wax has been removed, coat both halves of waxed area with a Cyanoacrylate or Diamond D Separator allow to dry and then close flask and tighten bolts.

10. Next, turn on the FlexPress™ and press the ON/STBY button. Flexpress may be turned on earlier before boil out stage to save time. For small tube, Program Flexpress to reach a heating temperature of

542°F with a melting time of 16 minutes and a hold time of 2 minutes and a recommended PSI of 105. For larger tube use a melting time of 19 minutes and a hold time of 2 minutes.

Important: Cut an “X” into the top of the ClearMet tube for easier injection. You may use a small disc to do this. After cutting an X into tube, Place small tube in toaster oven 175-200 degrees for 20 minutes. Place larger tube in toaster oven for 30 minutes. This will eliminate any moisture, which would affect the clarity and injection capabilities of the material.

11. When the FlexPress™ has reached the preset temperature, a buzzer will sound. At this time, insert the desired sized ClearMet cartridge into the heating chamber with the crimp end down. A light coat of release compound on the cartridge will help eliminate the cartridge getting stuck in the heating chamber. Press the start button on the Flexpress. This will start the injection cycle.

12. Place the flask into the FlexPress™, centering the sprue hole of the flask over the center of the heating chamber. Lock the flask in place using the flask holder and bolts.

13. When ClearMet has reached the required melting temperature, the FlexPress™ will automatically inject the ClearMet cartridge. The piston will stay in the up position according to the hold time. Once the piston retracts, it is safe to unscrew the hold down plate and remove the flask.

14. To remove the flask, loosen the hold down knobs and remove the flask holder. Using heat resistant gloves, grasp the flask and press the EXTEND/RETRACT button (lower right button). The will extend the piston and push the crushed cartridge and flask out for easy removal.

15. After removing the flask, press the EXTEND/RETRACT button again to retract the piston.

16. Cool flask for 10-15 minutes under cool water and the de-vest appliance. After the appliance is divested, it can be placed in an ultrasonic cleaner to eliminate any extra stone that is adhering to the appliance. The Cyanoacrylate or Diamond D Separator creates a smooth, clean surface on the underside of the appliance.

17. To create holes in the saddle areas use a standard round bur. This will allow the acrylic to mechanically retain to the ClearMet framework.

18. Fit, finish, polish with pumice or Hatho Thermopol Kit and then set teeth on frame for try in or for processing.

19. Before packing partial denture, use a bonding agent on mechanical retention areas.

20. Make sure to extend acrylicl onto the clasp where inter proximal areas are to avoid a black or grey triangle appearance.

21. When case is finished and polished coat with a very thin layer of Palaseal®. Then light cure for three minutes only.

NOTE: To tighten clasp LIGHTLY warm a 3 prong plier and lightly crimp clasp for 5 seconds.

For additional product information: see keystoneind.com

Please Note: We accept and assume no responsibility for consequential, incidental, special, indirect or punitive damages from the use of this product. Any responsibility for loss or damage arising out of or resulting from the use of this product shall not exceed the purchase price paid for the product.

FlexPress™ and Palaseal® are not registered trademarks of Keystone Industries.

keystone [®] <small>INDUSTRIES</small>	ECI/REP	Keystone Europe LLC Batavenweg 7 5349BC Oss, NL	CE
52 West King Street Myerstown, PA 19426 800-333-3131 www.keystoneind.com		NATLAB-1015 1980249	Rev 01/2017

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