

+ OPTIMIZE

AUTOCLAVE CLEANER

AUTOCLAVE CLEANER

LIQUID CONCENTRATE

Steam autoclave cleaner with enhanced cleaning power. Contains no phosphates or mercury. Cleans grease, grime, and particles from all Sterilizers. Powerful anionic and non-ionic surfactants.




Keep out of reach of children.

Net Contents:
32 fl oz (946 ml)

Makes 2.11 gallons (8 L)
1:16 dilution ratio

REF 00146



Optimize Autoclave Cleaner is a mildly alkaline solution used to remove grease, grime and scale deposits from autoclaves. Use as part of a routine maintenance program will keep autoclaves operating at peak performance.

Directions:

1. Dilute 4 fl oz (118.29 ml) of Optimize Autoclave Cleaner with 64 fl oz (1.89 L) of distilled water. For extremely dirty autoclaves, dilute 6 fl oz (177.44 ml) of Optimize Autoclave cleaner with 64 fl oz (1.89 L) of distilled water.
2. Drain reservoir and fill with prepared cleaning solution from Step #1.
3. Run one cycle without instruments in normal manner. Do not sterilize instruments while cleaning autoclave.
4. Drain cleaning solution from autoclave and reservoir and rinse until water runs clear. Fill reservoir again with distilled water and run one cycle.
5. Drain rinse solution from reservoir, remove tray rest and wipe out inside of chamber and tray rest.
6. Add distilled water to reservoir. Autoclave is now ready for use.

Warning

Causes serious eye irritation. Wear protective gloves/clothing/eye protection. Refer to product SDS for additional health, safety and disposal information.

Contains alkaline ingredients. Avoid contact with skin, eyes and clothing. In case of contact, flush with water. If irritation continues, obtain medical attention.

Do not take internally. Contains no phosphates or mercury.



Keystone Industries
52 West King St.
Myerstown, PA 17067 USA
800.333.3131
keystoneindustries.com

Made in USA

keystoneindustries.com/symbolglossary

International Sales

Keystone Industries GmbH,
Stockholzstr. 11, 78224 Singen, Germany
☎ +49 77 31 91 21 01

NATLAB-1359 4926639 02/2021