

+OPTIMIZE

ULTRASONIC SOLUTION

GENERAL PURPOSE ULTRASONIC SOLUTION NON-AMMONIATED POWDER CONCENTRATE




Keep out of reach of children.

Net Contents:
4.23 oz (120 g)
Makes 1 gallon (3.78 L)

REF 00158

CAUTION: US Federal law restricts this device to sale by, or on the order of a dental professional.



Advanced formulation combines nonionic surfactants, powerful chelators and corrosion inhibitors that enhance the cleaning action produced by ultrasonic cleaners while protecting and extending the life of delicate instruments. Effectively removes dried blood, saliva, pumice, rouge, buffing compounds, food debris, foreign matters and soils of all types from instruments and appliances when used in conjunction with an ultrasonic cleaner. May be used directly in the ultrasonic tank or as a sonic conductor when using Optimize Tartar & Stain Remover.

10:1 Concentrate provides for easier storage.

For use prior to final sterilization.

Directions:

1. Add warm water to shoulder of this container, replace cap, shake.
2. Total immersion time for items not to exceed 30 minutes.
3. Rinse items thoroughly with water prior to sterilization.
4. Prepare fresh solution daily or when solution becomes cloudy.



Caution: Alkaline solution can cause irritation upon contact with skin. Immediately flush eyes with water upon contact.

Danger

Causes serious eye irritation. Causes skin irritation. May be harmful if swallowed. May cause respiratory irritation if inhaled. Wear protected gloves/protected clothing/eye protection. Refer to product SDS for additional health, safety and disposal information.



Keystone Industries
52 West King St.
Myerstown, PA 17067 USA
800.333.3131
keystoneindustries.com

International Sales

Keystone Industries GmbH,
Stockholzstr. 11, 78224 Singen, Germany
+49 77 31 91 21 01

Made in USA

keystoneindustries.com/symbolglossary

EPR-40 NATLAB-1349 4926626 02/2021