

+OPTIMIZE

REFRESH

REFRESH MULTI-ENZYME GENERAL PURPOSE ULTRASONIC LIQUID CONCENTRATE

Directions:

For manual or ultrasonic cleaning, dilute 1 fl oz (29.57 ml) **Refresh** liquid concentrate with 1 gallon (128 fl oz, 3.78 L) of warm water. Cleaning activity increases with higher temperatures. Soak instruments for a minimum of 5 minutes. Soak instruments longer if heavily soiled or there is dried organic matter. Do not exceed 130°F (55°C). After cleaning, remove instruments from solution and thoroughly rinse with plain water. Appropriate sterilization and/or disinfection of instruments should follow cleaning process. Discard solution daily or sooner depending on soil level of instruments being cleaned.

Store in a cool place below 88°F (31°C).

Caution

Do not breathe mist, vapors or spray. Wear protective gloves. Wear protective clothing. Wear eye and face protection. In case of inadequate ventilation, wear respiratory protection.

Danger

Causes skin irritation. Causes eye irritation. May cause allergy or asthma symptoms or breathing difficulties if inhaled. May cause respiratory irritation. Refer to product SDS for additional health, safety and disposal information.



keystone®
industries

Keystone Industries
52 West King St.
Myerstown, PA 17067 USA
800.333.3131
keystoneindustries.com

International Sales

Keystone Industries GmbH,
Stockholzstr. 11, 78224 Singen, Germany
☎ +49 77 31 91 21 01



Made in USA

keystoneindustries.com/symbolglossary

NATLAB-1353 4926656 05/2021