# Bosworth Light Liner®

## Light Cure Chairside Reline Composite

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**Refer to product SDS for additional health, safety and disposal information.**

**Warning:** This product can expose you to chemicals including, N,N-dimethyl-p-toluidine and titanium dioxide, which are known to the State of California to cause cancer and dibutyl phthalate, which is known to the State of California to cause birth defects or other reproductive harm. This product can also expose you to cadmium, which is known to the State of California to cause cancer, and birth defects or other reproductive harm. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

### Liquid

#### Danger

Harmful if inhaled (Applies to Hard Reline Liquid Only). Causes serious eye irritation. Causes skin irritation. May cause an allergic skin reaction. May damage the unborn child. Suspected of damaging fertility. May cause respiratory irritation (Applies to Hard Reline Liquid Only).

### POWDER

#### Warning

May cause an allergic skin reaction. May cause cancer.


### AIR BARRIER

#### Warning

Harmful if swallowed. Causes eye irritation.

Keystone Industries  
52 West King St.  
Myerstown, PA 17067 USA  
800.333.3131  
[keystoneindustries.com](http://keystoneindustries.com)

EC REP  
Keystone Europe LLC  
Batavenweg 7  
5349BC Oss Netherlands

Liquid Flashpoint 230°F (110°C)

[keystoneindustries.com/symbolglossary](http://keystoneindustries.com/symbolglossary)

**Made in USA**

PR098 BosLab-309 08/2020


# Bosworth LightLiner Instructions for Use

Bosworth LightLiner is a unique light-cure chairside reline material that is available in a soft or hard non-methyl methacrylate formula.


## PREPARATION

Remove denture from patient's mouth and grind 1-2mm from the entire palatal surface with a bur to expose fresh acrylic base. Coat patient's teeth and labial surface with a silicone lubricant such as Bosworth Masque to prevent the composite from sticking.


## APPLICATION


1. Add 1 vial of liquid (4g) to plastic cup of powder (6g). Spatulate for 30 seconds to obtain uniform mix. Mix slowly to reduce bubbles.

2. Coat denture with a thin layer of LightLiner mix to serve as a bonding agent.<sup>1</sup>

3. Set aside remainder of mix until it reaches a honey-like consistency (2-3 minutes).

4. Spatulate mix into denture and place in patient's mouth. Wipe off excess. Hold in patient's mouth until composite reaches a rubber-like consistency (3-5 minutes).

## CURING

### Soft Formula

1. Remove denture from patient's mouth, rinse and dry.
2. Place in a light cure unit with a 100-watt halogen for 5 minutes.<sup>2</sup>
3. Remove denture and coat with Air Barrier.
4. Cure for an additional 5 minutes.
5. Wash denture to remove Air Barrier coating.
6. Trim and polish as necessary.

### Hard Formula

1. Remove denture from patient's mouth, rinse and dry.
2. Coat denture with Air Barrier.
3. Place in a light cure unit with a 100-watt halogen for 5 minutes.<sup>2</sup>
4. Wash denture to remove Air Barrier coating.
5. Trim and polish as necessary.

## NOTES FOR USING LIGHTLINER

1. If air bubbles surface, puncture immediately with a spatula to eliminate any porosity.
2. Times listed above are applicable only if curing with a 100 watt oven light. If curing with a 250 watt oven light, time will be approximately 1 1/2 minutes. May be cured with a hand-held unit with a wavelength range greater than 450 nm (approximately 10-15 minutes). A more powerful light source will reduce the time of cure and a less powerful light source will increase the time of cure.

CAUTION: Certain individuals may be sensitive to products containing acrylic.

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