Diamond D® Heat Cure Processing Instructions:

Keystone Industries' Diamond D Heat Cure Denture Acrylic is intended to be used in the fabrication of removable and permanently fixed partial and full dentures. It is certified as Class IIa in Europe and is approved for use with implants.

1. Flask

Wax up, flask, and boil out in the usual manner. Make sure all surfaces are wax free before proceeding. For mold washout, we recommend Diamond D Ultra Wash. When steam has stopped rising from the cast, apply Diamond D Ultra-Sep Tinfoil Substitute to all mold surfaces. Avoid coating teeth with this material since this can result in poor adhesion between teeth and the denture acrylic. If teeth are coated with tin foil substitute, remove it from the teeth with the Diamond D Sep Brush II and a little monomer.

2. Mixing

Slosh/roll powder container before each use to insure proper fiber dispersion. **Powder/Liquid Ratio by weight:** 10 ml monomer to 21g powder. **Powder/Liquid Ratio by volume:** 10 ml monomer to 30 cc powder.

For optimal results follow the above ratios and process at ambient room temperature. Measure liquid and pour into a glass mixing jar with lid. Shake the powder to evenly distribute the fibers. Measure powder and add enough to the liquid while stirring until the liquid is completely absorbed. Tap the jar on the bench top to bring the remaining monomer to the surface and add any remaining powder. Stir with a stainless steel spatula until all the powder is moistened. Cover and let the acrylic rest until it pulls away from the sides of the jar and does not stick to the stainless steel spatula. Estimated time for Diamond D Heat Cure Acrylic to reach packing stage is 7 minutes to 14 minutes depending on ambient air temperature.

3. Packing

DO NOT LET ACRYLIC REACH "SNAP STAGE". Acrylic that is packed too dry can result in brittleness, porosity and shade discoloration. Estimated working time after the material has reached packing stage is 5 to 7 minutes. The flask should be at room temperature 60° - 80°F (15° - 27°C) prior to packing. Remove the acrylic from the mixing jar. Trial pack denture 2-3 times to insure minimized porosity and maximum acrylic density. Use sufficient pressure with each trial pack and remove excess acrylic using polyethylene sheets as a separator between each trial pack. Remove the flask and close final without the polyethylene sheet. For best result allow clamped flask to sit on bench or in room temp water bath for 20-30 minutes to allow residual internal exothermic heat to subside, then begin curing process.

4. Processing

Water temperature is critical in the processing of all denture acrylics. Polymerization generates high temperatures inside the flasks. If the flask is put immediately into boiling water the liquid resin may boil and cause porosity in the denture base. The curing cycles below insure porosity free castings. For thicker cases bench set in flask compress for 30 minutes before introducing into 165°F (73°C) water bath and process through the normal curing cycle.

It has been noted in many clinical studies that longer cure times produce less residual monomer and a terminal boiling treatment for at least 1 hour is known to achieve maximum monomer conversion. An additional hot water bath for 1 hour at 50°C prior to insertion of the device into the oral cavity is also recommended for minimizing the possible risk of allergic reactions to the denture acrylic.

Long Cure - Immerse flask into room temperature water, rise to 165°F (73°C) and maintain that temperature for 8 - 10 hours.

Heat polymerization/long cure Option 2

After 20 minute resting period place packed or injected and clamped flask into room temperature water bath and raise temperature to 165°F (73°C). Hold at temperature overnight or 8 - 10 hours. If desired an optional boil for final 60 minutes of long cycle is possible. Remove from bath and allow flask to cool on the counter for 20 minutes before final immersion in lukewarm water for final cool. Deflask and finish.

Short Cure - Immerse flask into room temperature water, rise to 165°F (73°C) and maintain that temperature for 90 minutes. Then raise temperature to a boil 212°F (100°C) and hold for 30 minutes. Allow to cool and then deflask and finish. **Note: To reduce any residual monomer which may result in the risk of patient sensitivity, you may leave flask in boiling water for an additional 30 minutes.**

Heat polymerization/short cure Option 2

After 20 minute resting period place packed or injected and clamped flask into room temperature water bath and raise temperature to 165°F (73°C). Hold for 90 minutes, and then raise temperature to a boil 212°F (100°C) and hold for 30 minutes. Remove from bath and allow flask to cool on the counter for 20 minutes before final immersion in lukewarm water for final cool. Deflask and finish. **Note: To reduce any residual monomer which may result in the risk of patient sensitivity, you may leave flask in boiling water for an additional 30 minutes.**

5. Deflasking

When the flask is cool enough to handle, place it under warm water for 20 minutes or bench cool for 45 minutes (recommended) then deflask in the usual manner. In the case of undercuts, do not pry the denture with a knife or distortion of the denture could result. Section the cast instead to insure a superior fit. Keep unfinished denture in a moist environment until it can be finished.

6. Finishing and prosthesis storage

After post process equilibration, prosthesis is removed from cast and trimmed with conventional crosscut burs. Prior to pumicing, miniature Scotchbrite wheels can provide a pre-pumice surface allowing for ideal final polishing. It is always recommended to store prosthesis in water during the course of finishing and until it is inserted in mouth.

7. Product Storage

Liquid should be stored in a temperature constant area not to exceed 80°F (27°C). Avoid storage in direct sunlight.

*Always practice Good Manufacturing Practice protocols regarding the non-mixing of powders or liquids of one manufacturer with another.

NOTE: When tested according to ADA/ANSI Spec. #12, 8.7, the maximum residual monomer content of Diamond D heat cure material, when processed as outlined above, will be less than 2%.

Precautions:

Product contains methyl methacrylate monomer which, although rare, may cause an allergic reaction in individuals sensitive to acrylic containing products. Flash point is 50°F (10°C). Keep away from heat, sparks and open flame. Avoid inhalation or ingestion. Use in a well-ventilated area. Keep container closed when not in use. Avoid prolonged or repeated contact with skin. Wash skin with soap and water after contact. Keep out of reach of children. Use only non-alcoholic separators with this acrylic.

Keystone Industries
52 West King Street
Myerstown, PA 17067 USA
800-333-3131
keystoneindustries.com
sales@keystoneind.com

CE REP
Keystone Europe BV
Bataveweg 7
5349 BC Oss, NL

Keystone Industries GmbH,
Stockholzstr. 11, 78224 Singen, Germany
+49 77 31 91 21 01

Report any serious incident occurring with this device to the manufacturer and applicable Competent Authority of the member state in which the user/patient is established.

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