Cemento FLECK'S®

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FLECK'S® CEMENT

Oxyphosphate of Zinc - Fine Particle Size

Complies with ISO 9917-1 for dental water-based cements.

Indications for Use: Fleck's Cement is indicated for use as a luting agent to seal dental appliances (e.g., crowns, bridges, facings, and orthodontic brackets) to teeth or to other appliances, as a base or liner for restorative materials, and as a tempo-rary restorative material in deciduous teeth.

Warnings/Precaution: Liquid can cause severe skin burns and eye damage. Powder causes serious eye irritation. Although rare, hypersensitivity reactions have been noted with the use of other zinc oxide products; therefore, the possibility exists for developing a local sensitivity or allergic response to Fleck's Cement.

Storage: Keep containers tightly closed and protected from direct sunlight when not in use and store in a dry area at normal room temperature away from volatile material.

Disposal: Waste is harmful to aquatic life with long lasting effects and should not enter sewage system. Dispose of unused materials in accordance with federal, state, and local regulations.

Mixing Utensils: One clean, dry glass mixing slab at a temperature of 18-24°C (65-75°F), a corrosion resistant spatula. The slab and spatula must be free from dust, set cement, moisture, etc.

Mixing Technique for General Use: Cement preparation should be done at ambi-ent room temperature (65-75°F). 0.8 grams of powder (2 fillings to the shoulder of the powder cap dome) are placed on the slab and formed into a flat rectangular shape with a spatula. The powder is divided first into quarters, one of these quar-ters is divided into 2-eighths, and one of these eighths is divided into 2-sixteenths, so that there are 6 portions; 3-quarters, 1-eighth and 2-sixteenths. 0.3 ml of liquid (12 drops) is placed on the slab; add 1-sixteenth portion of powder to the liquid and mix with a linear motion for 15 seconds. Add the 2nd sixteenth portion, mix for 15 seconds, then the 1-eighth portion, mix for 15 seconds, followed by 1-quarter portion, mix for 20 seconds, the 2nd quarter portion, mix for 20 seconds and the last quarter, mix for 35 seconds, completing the mix in two minutes. Each portion is to be thoroughly incorporated, with no unmixed particles of powder or liquid remaining around the edges.

A linear, not rotary, motion of the spatula is to be used with the edge of the spatula sweeping approximately half of the mixing area of the slab on each stroke. A mix made in this manner will flow reluctantly from the spatula and be so completely plastic as to allow the accurate placement of an inlay or crown no matter how precisely fitted. Maximum working time is 3-1/2 minutes. Net setting time is 2-8 minutes.

Finishing: Excess material may be removed at once or when the cement has com-pletely hardened using a hand held instrument.

Mixing Technique

Liquid 0.3ml powder 0.8g	20 secs.	20 secs.	20 secs.
	1/4	1/4	1/4
Mixing time of each portion in seconds. Total mixing time is 120 secs.	15 secs.	15 secs.	15 secs.
	1/16	1/16	1/8

Thick Mix: Obtain a thick mix by using a larger portion of powder with the same procedure.

Setting Time Control with the One Liquid: SLOW SETTING: Spatulate a pin-head of powder into the liquid required for your mix. Let this stand 2 or 3 minutes. proceed to make your mix in the regular manner, spatulating as described above. **QUICK SETTING:** Reduce mixing time by reducing the mixing time of each portion by half (see diagram), total mixing time, 1 minute.

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Keystone Industries
 52 West King Street, Myerstown, PA 17067 USA

(800) 333-3131 | keystoneindustries.com

To order outside of the US call +49 (0) 77 31 91 21 01

EC RE P

 Keystone Europe BV
 Batavenweg 7

 5349 BC Oss, NL

CE
1 639
 MD