

Section 7	Handling and storage
Handling:	Wear safety glasses with side shields (recom: ANSI Stand Z87.1)
Measures against fire and explosion:	Use exhaust for polishing dust
Storage	In dry surrounding
Section 8	Exposure controls/personal protection
General:	Use general rules for polishing with rotating instruments
Personal protection:	Wear eye protection and gloves when necessary
Respiratory protection:	Avoid breathing of dust created by cutting, sanding or grinding. Select one of the following NIOSH approved respirators based on airborne concentration of contaminants and in accordance with OSHA regulations: half-Mask dust respirator, full-face dust respirator.
Eye protection:	Wear safety glasses
Section 9	Physical and Chemical Properties
Melting point:	> 35 °C (104°F)
Boiling point:	> 130 °C (300°F)
Relative density:	
Solubility:	Insoluble in water
Colour	Blue
Appearance	Pasty
Section 10	Stability and Reactivity
Conditions to avoid:	None
Materials to avoid:	None known
Hazardous Polymerization:	Will not occur.
Hazardous decomposition products:	Incomplete combustion will produce smoke, CO and CO ₂ .
Section 11	Toxicological information
Routes of entry:	
Inhalation:	Due to the low vapor pressure inhalation is unlikely under normal temperatures. Inhalation of decomposition products like vapor or smoke may cause adverse health effects
Eyes:	This product is abrasive. Prolonged contact may cause temporary skin and eye irritation.
Section 12	Ecological Information
	No effect
Section 13	Disposal Considerations
	None
Section 14	Transport Information
	No special care necessary
Section 15	Regulatory Information
	None
Section 16	Other Information
<p>The purpose of this safety sheet is to describe the product in terms of health and safety and not as a product specification, guaranteeing their properties.</p> <p>The information on this Data Sheet is accurate to the best of our knowledge as to the proper use and handling of this product under normal conditions. Any use of the product which is not in conformance with this Data Sheet or which involves using the product in combination with other products or any other process is the responsibility of the user.</p> <p>HATHO makes no warranty, either express or implied, with respect to the completeness or continuing</p>	

accuracy of the information contained herein and disclaims all liability for reliance thereon.