



Eugenol USP

About Eugenol USP

Eugenol is a non-irritant liquid to be mixed with Zinc Oxide powder in order to fabricate a low strength Zinc Oxide Eugenol (ZOE) temporary cement.

Features & Benefits

- Intended to serve as a temporary tooth filling or as a base cement to affix a temporary tooth filling.
- Is used to affix dental devices such as crowns or bridges.
- Can be applied to a tooth to protect the tooth pulp.
- Made in the USA



PRODUCT

QUANTITY

ITEM NUMBER

Eugenol

1 oz (30 mL) bottle

#06-02077

keystone[®]
industries

dental.keystoneindustries.com

US: 480 South Democrat Road Gibbstown NJ, 08027
EU: Keystone Industries GmbH Stockholzstr. 11 78224 Singen, Germany

US: 1.800.333.3131
EU: +49 77 31 91 21 01