



Hatho Miniature Scotch Brite™

Hatho Miniature Scotch Brite Brushes are used to polish without ruining delicate anatomical surfaces on dental appliances.

About Scotch Brite™

By using low speeds – around 5,000 rpm – and light pressure the acrylic should not overheat due to friction. These brushes come in coarse, medium and fine to let you control the type of polishing finish you need.



Features & Benefits

- Available in Coarse, Medium, Fine
- Available in 12/Pkg

PRODUCT	QUANTITY	ITEM NUMBER
Coarse	12/Pkg	1670016
Medium	12/Pkg	1670092
Fine	12/Pkg	1670091

keystone[®]
industries

dental.keystoneindustries.com

US: 480 South Democrat Road Gibbstown NJ, 08027
EU: Keystone Industries GmbH Stockholzstr. 11 78224 Singen, Germany

US: 1.800.333.3131
EU: +49 77 31 91 21 01