



Hemox-A

Buffered Hemostatic Solution

About Hemox-A

Helps with crown and bridge cases for use with the retraction cord to control bleeding and displace gingival tissue for better impressions. A plain gingival retraction cord can be dipped in Hemox-A Solution and then placed into the gingival sulcus for its hemostatic action.

Features & Benefits

- 25% buffered Aluminium Chloride solution
- Contains no epinephrine
- Great for gingival retraction and hemostasis
- Made in the USA



PRODUCT	QUANTITY	ITEM NUMBER
Hemox-A	1 oz (30 mL) bottle	#25-03677



dental.keystoneindustries.com

US: 480 South Democrat Road Gibbstown NJ, 08027
 EU: Keystone Industries GmbH Stockholzstr. 11 78224 Singen, Germany

US: 1.800.333.3131
 EU: +49 77 31 91 21 01