

# Hinged EZ Fluoride Trays

## Fluoride Trays

These Hinged EZ Fluoride Trays are an excellent way to dispense fluoride treatments with upper and lower arches attached together.

The tray is made of closed celled, 4 lb. density foam for good rigidity. With excellent anatomical depth for gingival line, the tray promotes full coverage of the fluoride applications.



### FEATURES AND BENEFITS

- Sizes: Large, Medium, Small/Medium, Small.
- Sizes are color coded with four different colors representing four sizes.
- Placement handle for easy insertion.
- Available in 100 to 50 in a bag.

### ORDERING INFORMATION

Large – 100/bag	#9524860
Large Dual – 100/bag	#9525050
Medium – 100/bag	#9525200
Medium Dual – 100/bag	#9525390
Small/Medium Dual – 100/bag	#9525890
Small – 100/bag	#9525520

