

Edentulous Tray Aways®


Single-use, perforated edentulous impression trays that can be easily drilled for implant impressions.

About Edentulous Tray Aways®

These unique and economical trays allow for maximum material retention and accuracy without the use of an adhesive. These high impact styrene trays provide an excellent impression with outstanding details. Edentulous Tray Aways are anatomically designed to be patient friendly and to allow for a perfect fit in just one visit.









Features & Benefits

- Made from high impact styrene
- Perforated for maximum material retention
- Adhesive optional
- Easily adapted for implant impressions
- Ergonomic with built in handle
- Economical and disposable
- Anatomically designed
- No overflow of impression material
- Bur trimmable
- Non-toxic, no patient reaction



keystone[®]
Industries

dental.keystoneindustries.com

	ITEM NUMBER	SIZE	DESCRIPTION	QUANTITY*	DIMENSIONS
	#0921600	S	#1 Upper	12/bag	(A) 3.06" (B) 2.41" (C) 0.88"
	#0921887	L	#2 Perforated Lower	12/bag	(A) 3.02" (B) 2.33" (C) 0.86"
	#0921885	M	#3 Perforated Upper	12/bag	(A) 2.91" (B) 2.21" (C) 0.83"
	#0921888	M	#4 Perforated Lower	12/bag	(A) 2.93" (B) 2.18" (C) 0.79"
	#0921886	S	#5 Perforated Upper	12/bag	(A) 2.66" (B) 2.00" (C) 0.85"
	#0921889	S	#6 Perforated Lower	12/bag	(A) 2.78" (B) 2.06" (C) 0.77"
	#0921890	PARTIAL	#7s Solid UR, LL	12/bag	(A) 2.24" (B) 2.29" (C) 0.79"
	#0921891	PARTIAL	#8s Solid UL, LR	12/bag	(A) 2.15" (B) 2.33" (C) 0.86"