

Supergel®



Supergel® is a premium blended alginate material that mixes readily to a smooth, creamy consistency.

About Supergel®

Supergel's slow dehydration and low surface tension provide consistent, accurate, and bubble-free detailed reproductions every time. Pleasing sassafras flavor, controlled flow, heavy body, and gag reflex reduction provide greater patient acceptance.

Instructions

Fluff powder by rotating or shaking bag before using. Use 72°F (23°C) water. If necessary to delay the pouring of the model, wrap a wet towel around the impression and store in a plastic bag or sealed container.



Features & Benefits

- Smooth, creamy consistency
- Bubble-free reproductions
- Sassafras flavor
- Gag reflex reduction

PRODUCT	SIZE	ITEM NUMBER
Regular Set, Can	1lb (454g)	0921825
Fast Set, Can	1lb (454g)	0921826
Regular Set, Pouch	1lb (454g)	0921818
Fast Set, Pouch	1lb (454g)	0921819



dental.keystoneindustries.com