

No. 1
SELLING U.S.
TRAY



Tray Aways®



The Tray Aways® Disposable Impression Trays are the original plastic tray!

About Tray Aways®

Tray Aways are the disposable impression trays of choice suited for use with alginate or other types of impression materials where detailed reproduction and dimensionally-stable results are vital. Made from rigid, high-impact FDA approved styrene, Tray Aways assist in generating a consistent, unsurpassed impression free of distortion or bite through. Smooth, rounded edges, interior retention slots, and vertical perforation enhance material retention. The best value on the market.













Features & Benefits

- Made from rigid, high-impact styrene
- Produces dimensionally-stable results
- Creates impression free of distortion & bite through
- Also used with implants



keystone[®]
i n d u s t r i e s

dental.keystoneindustries.com

	ITEM NUMBER	SIZE	DESCRIPTION	QUANTITY*	DIMENSIONS
	#0921884	L	#1 Perforated Upper	12/bag	(A) 3.06" (B) 2.41" (C) 0.88"
	#0921887	L	#2 Perforated Lower	12/bag	(A) 3.02" (B) 2.33" (C) 0.86"
	#0921885	M	#3 Perforated Upper	12/bag	(A) 2.91" (B) 2.21" (C) 0.83"
	#0921888	M	#4 Perforated Lower	12/bag	(A) 2.93" (B) 2.18" (C) 0.79"
	#0921886	S	#5 Perforated Upper	12/bag	(A) 2.66" (B) 2.00" (C) 0.85"
	#0921889	S	#6 Perforated Lower	12/bag	(A) 2.78" (B) 2.06" (C) 0.77"
	#0921890	PARTIAL	#7s Solid UR, LL	12/bag	(A) 2.24" (B) 2.29" (C) 0.79"
	#0921891	PARTIAL	#8s Solid UL, LR	12/bag	(A) 2.15" (B) 2.33" (C) 0.86"
	#0921892	ANTERIOR	#9s Solid Lower	12/bag	(A) 2.56" (B) 1.50" (C) 0.75"
	#0921893	ANTERIOR	#10s Solid Upper	12/bag	(A) 2.52" (B) 1.47" (C) 0.76"
	#0921894	PARTIAL	#7 Perforated UR, LL	12/bag	(A) 2.24" (B) 2.29" (C) 0.79"
	#0921895	PARTIAL	#8 Perforated UL, LR	12/bag	(A) 2.15" (B) 2.33" (C) 0.86"
	#0921896	ANTERIOR	#9 Perforated Lower	12/bag	(A) 2.56" (B) 1.50" (C) 0.75"
	#0921897	ANTERIOR	#10 Perforated Upper	12/bag	(A) 2.52" (B) 1.47" (C) 0.76"

*Can be bagged in quantities other than 12.