

Acid Etch®

Acid Etch® is 35% high viscosity phosphoric acid. Our Acid Etch® is silica and glycerin-free, providing a deeper etch pattern for easier application of bonding material. For increased patient comfort and easy clean-up, Acid Etch® is highly viscous, allowing for more accurate placement on enamel or dentin. Material will not slump, slide or run. Acid Etch® is clear aqua blue in color, odorless and tasteless. Allows flowability and clarity while providing complete coverage. Ideal for everyday restorative procedures within every office.



Jumbo Syringe Kit



Standard Kit


ORDERING INFORMATION

ACID ETCHANT

0921303 Jumbo (60 gm) Syringe - 35%



0921346 Standard Kit - 35%



0921345 1.2 ml Kit



RECOMMENDED FOR USE WITH ACID ETCH

GORILLA SEAL™

Ultra-clear light-cure Sealant

AEGIS® PIT & FISSURE

Pit & Fissure Sealant with ACP

FEATURES AND BENEFITS

- 1 Clear aqua blue in color
- 2 Silica-free & glycerin-free
- 3 Deeper etch pattern for effortless application of bonding material
- 4 Highly viscous material enables more accurate placement
- 5 Odorless & tasteless for ultimate patient comfort
- 6 Easily removed with a high speed evacuator