

# Edentulous Tray Aways®

Edentulous Tray Aways® are single-use, perforated edentulous impression trays that can be easily drilled for implant impressions. These unique and economical trays allow for maximum material retention and accuracy without the use of an adhesive. These high impact styrene trays provide an excellent impression with outstanding details. Edentulous Tray Aways are anatomically designed to be patient friendly and to allow for a perfect fit in just one visit. Offered in three convenient sizes: small, medium, and large.



## ORDERING INFORMATION

### EDENTULOUS TRAY AWAYS

0921600	#1 Upper Small (Bag of 12)
0921599	#2 Lower Small (Bag of 12)
0921700	#3 Upper Medium (Bag of 12)
0921701	#4 Lower Medium (Bag of 12)
0921702	#5 Upper Large (Bag of 12)
0921703	#6 Lower Large (Bag of 12)
0921706	#7 Upper X-Large (Bag of 12)
0921707	#8 Lower X-Large (Bag of 12)
0921704	Variety Bag (2 of each size above)

CAN BE BAGGED IN QUANTITIES OTHER THAN 12

### RECOMMENDED FOR USE WITH AIR POLISHER

**BOSWORTH SUPERGEL®**  
Alginate Impression Material

**BOSWORTH SUPERSIL™**  
Polyvinylsiloxane Impression Material

**BOSWORTH TAC®**  
Tray Adhesive Compound

## FEATURES AND BENEFITS

- 1 Made from high impact styrene
- 2 Perforated for maximum material retention
- 3 Adhesive optional
- 4 Easily adapted for implant impressions
- 5 Ergonomic with built in handle
- 6 Economical and disposable
- 7 Anatomically designed
- 8 No overflow of impression material
- 9 Bur trimmable
- 10 Non-toxic, no patient reaction