


Esteem®

Esteem® is a vita shade (A-2), 100% radiopaque, self-cure composite restorative material. The two-paste system, which includes a universal paste and catalyst paste, forms a durable restoration. The hardened surface can be polished to a natural appearance, almost identical to original tooth anatomy. Esteem® is stain resistant and has physical properties that produces exceptional holding strength. It matches the full range of tooth structure with a single translucent shade that blends into the surrounding dentition. It mixes easily, sets quickly, and finishes beautifully. Can be used for post and core build-up as well.

ORDERING INFORMATION

ESTEEM®

0921560 Standard Kit



0921564 Double-Ended Spatulas (100)



FEATURES AND BENEFITS

- 1 Two paste system
- 2 Radiopaque
- 3 Natural, luster look
- 4 Available in orange, black, red, blue
- 5 Durable and translucent
- 6 Contours and polishes easily

Gorilla Seal™

Gorilla Seal™ is a 100% crystal-clear light-cure surface sealant that gives a gloss coat to composite restorations and temporaries. May reduce staining and sensitivity with bleaching, cervical erosion, and veneers. Gorilla Seal™ can be used to seal over I.D. or decorative designs.


ORDERING INFORMATION

GORILLA SEAL™

0921050 Standard Kit (bottle)



0921055 Standard Kit (syringe)



Standard Kit (bottle)



Standard Kit (syringe)

FEATURES AND BENEFITS

- 1 Crystal clear
- 2 Cures to a gloss finish allowing for a reduced polishing
- 3 May reduce sensitivity when used in cases such as tooth beaching, cervical erosion and veneers

Hydrox[®]

Hydrox[®] radiopaque calcium hydroxide paste is a rigid, self-setting cavity liner and pulp capping agent. Hydrox[®] eliminates all chemical and thermal sensitivity by sealing dentinal tubules and protecting the pulp. Hydrox[®] serves as a protective barrier that insulates the dentin and pulp from acid-containing cements and restorative materials. Works well under glass ionomers and will not affect the amalgam or the polymerization of composites. The superior flowability of Hydrox[®] permits uniform coverage and controlled placement into the preparations.



ORDERING INFORMATION

HYDROX[®]

0921071 Standard Kit, Dentin

0921072 Standard Kit, Ivory



0921071 Standard Kit, Dentin

0921072 Standard Kit, Ivory



FEATURES AND BENEFITS

- 1 Easy to place & mix
- 2 Great temporary cement
- 3 Eliminates all sensitivity
- 4 Cavity Liner and Pulp capping agent
- 5 Provides protective barrier
- 6 Self-curing