

Hatho ThermoPol Set

Less work for better results is what the new Hatho ThermoPol Set is built upon. It's revolutionary design and working technique will brighten not just acrylics in faster-than-ever times, but also thermoplastics such as Valplast® and TCS®.

A combination of the Hatho Multi Layer Brush and Polistar Lintgrey Compound is comparable to no two other products for the pre-polishing process. The Multi Layer Brush is equipped with specially selected cloth or fiber in a design that retains the Compound on the brush, leaving you with an evenly polished thermoplastic after three minutes. Thanks to the high-concentration polishing abrasives found in Polistar Cream and Pro Pink that are used on the COSIMA and MIRA buffs, the sparkle and shine on the appliance will be brighter than ever before.



Hatho ThermoPol Set



Multi Layer Brush



Cosima Soft Buff



Mira Polishing Buff



Polistar Cream



Polistar Pro Pink



Polistar Lintgrey Compound

ORDERING INFORMATION

HATHO THERMOPOL SET

1670540	Hatho ThermoPol Set
1670541	Hatho Multi Layer Brush
1670542	Hatho Polistar Lintgrey Compound
1670543	Hatho Cosima Soft Buff
1670544	Hatho Polistar Cream
1670546	Hatho Polistar Pro Pink
1670547	Hatho AlphaCard Pocket Card
1670523	Hatho Mira Polishing Buff

FEATURES AND BENEFITS

- 1 Less work for better results
- 2 Will brighten thermoplastics such as Valplast® & TCS®
- 3 All-in-one kit is the only of its kind on the market
- 4 Quick three step process resulting in a happy patient