

# Duz-All®

**Duz-All is a PMMA methyl methacrylate all-purpose, self-curing acrylic repair material intended for any emergency.**

## About Duz-All®

Duz-All is intended for emergencies including: replacing broken teeth, repairing fractures, adding clasps, cold-cure full dentures, orthodontic retainers, rebasing and correcting full or partial dentures, adding a post-dam and appliances. It is perfect to use in a denture duplicator.

Sets in less than 20 minutes with no special equipment needed. Ideal for making or repairing jacket crowns.

## Features & Benefits

- Bulk packaging available in 5lb, 10lb, 25lb and custom



## Standard Kit ( 8 oz. Powder, 4 oz. liquid)

PRODUCT	ITEM NUMBER
Standard Kit, Pink	166260
Standard Kit, Clear	166260C
Standard Kit, Vein	166260V
Standard Kit, White	166260W

## Components & Bulk

PRODUCT	SIZE	ITEM NUMBER
Powder, Pink	8 oz. (227g)	166261
Powder, Clear	8 oz. (227g)	166261C
Powder, Vein	8 oz. (227g)	166261V
Powder, White	8 oz. (227g)	166261W
Liquid	4 oz. (118ml)	166262
Liquid Only	8 oz. (236ml)	166265
Liquid	32 oz. (946ml)	166269



[dental.keystoneindustries.com](http://dental.keystoneindustries.com)