



# Diamond D<sup>®</sup> Denture Acrylic

**Advanced natural-looking dentures requires artistry as well as technical expertise.**

## About Diamond D<sup>®</sup>

Since its introduction in 2003, Diamond D<sup>®</sup> has become one of the world's leading high-impact denture acrylic. Laboratories and dentists alike have come to rely on Diamond D's superior strength and natural fibered colors to deliver a premium, fracture-resistant denture product to patients.


## Features and Benefits

- Powder & liquid available in industry standard colors
- High-impact strength
- Exceptional flexural strength
- Exceeds expectation in patient comfort
- Consistent & forgiving
- Available in heat cure & self cure

IMAGE	COLOR	POUNDS	SKU
	ORIGINAL	1 LB	9576540
		5 LBS	9576530
		25 LBS	9576320
	LIGHT	1 LB	9576310
		5 LBS	9576360
		25 LBS	9576350
	LIGHT REDDISH PINK	1 LB	9576310
		5 LBS	9576360
		25 LBS	9576350
	CHROMA ESSENCE	1 LB	9576310
		5 LBS	9576360
		25 LBS	9576350
	DARK VEINED	1 LB	9576310
		5 LBS	9576360
		25 LBS	9576350

## Dental Acrylic Izod Impact Strength

Impact strength - a measure of toughness. The stress to fracture a notched specimen with a single blow. Expressed in foot-pounds of energy absorbed. Designed as "Charpy" or "Izod" impact strength depending on the testing machine used.



Source: Notched Izod Impact Strength Determinations per ASTM D256. Submitted in ASTM D638, Type 1 specimen form. Performed at Bodycote Polymer-Broutman Laboratory, a UL and IAPMO approved independent testing laboratory, A2LA and IAS accredited. BBI Test Reference #: 30002570 and 30000666.



[dental.keystoneindustries.com](http://dental.keystoneindustries.com)