

Fastray™

Fastray is a dimensionally stable, radiopaque, self-curing individual tray and acrylic base plate material available in regular (5 min.) or extra-fast (3 min.) set.

Features & Benefits

- Dimensionally stable, radiopaque, self-curing
- Available in regular (5 min.) or extra-fast (3 min.) set.
- No memory and when used with a separator will not stick to fingers or the model
- Available in 3 colors: Blue, Pink and White
- Recommended for use with Masque™



Bulk Kit (5 lb powder, 1 qt liquid)

PRODUCT	ITEM NUMBER
Bulk Kit, Regular Blue	0921386
Bulk Kit, Regular Pink	0921387
Bulk Kit, Regular White	0921388
Bulk Kit, X-Fast Blue	0921389
Bulk Kit, X-Fast Pink	0921390
Bulk Kit, X-Fast White	0921391

Components & Bulk

PRODUCT	SIZE	ITEM NUMBER
Powder, Blue	5 lb (2.2kg)	0921392
Powder, Pink	5 lb (2.2kg)	0921393
Powder, White	5 lb (2.2kg)	0921394
Liquid Regular	1 qt. (946ml)	0921395
Liquid X-Fast	1 qt. (946ml)	0921396

Standard Kit (1lb powder, 8 oz. liquid)

PRODUCT	ITEM NUMBER
Standard Kit, Regular Blue	0921375
Standard Kit, Regular Pink	0921376
Standard Kit, Regular White	0921377
Standard Kit, X-Fast Blue	0921378
Standard Kit, X-Fast Pink	0921379
Standard Kit, X-Fast White	0921380

Components & Bulk

PRODUCT	SIZE	ITEM NUMBER
Powder, Blue	1 lb (454g)	0921381
Powder, Pink	1 lb (454g)	0921382
Powder, White	1 lb (454g)	0921383
Liquid Regular	8 oz (236ml)	0921384
Liquid X-Fast	8 oz (236ml)	0921385
Powder, Blue	25 lb (11.34kg)	0921429B
Powder, Blue	100 lb (11.34kg)	0921428B



dental.keystoneindustries.com