

Hatho Geza Brushes

Hatho Geza's exclusive geared design and premium bristle materials offer the perfect miniature brush for deep, interproximal and occlusal polishing.









About Hatho Geza Brushes

Hatho Geza's exclusive geared design and premium bristle materials (Goat Hair and Chung King) offer the perfect miniature brush for deep, interproximal and occlusal polishing. It's a 3/32" shank and ideal speed is set at 15,000 RPM. Available in 12/Pkg.



Features & Benefits

- Ideal speed is set at 15,000 RPM
- A 3/32" shank
- Premium bristle materials (Goat Hair and Chung King)

PRODUCT PHOTO		QUANTITY	ITEM NUMBER	
		12/PKG	1670087	WHITE GOAT HAIR
		12/PKG	1670089	WHITE CHUNG KING
		12/PKG	1670088	BLACK CHUNG KING
		12/PKG	1670090	GREY CHUNG KING



dental.keystoneindustries.com