



5,000 rpm

Hatho Habras Discs

Habras Discs offer a unique, three-dimensional abrasive/flexible bristle design that provides a uniform finish.

About Hatho Habras Discs

Habras Discs provide the convenience of a pre-charged abrasive but do not generate the heat associated with rubber wheels. Embedded cubitron mineral abrasive is much harder than most metal. Removes oxides and polish without using compounds. Ideal for cleaning metal work and hard-to-reach areas. Ideal speed: 10,000 rpm. Discs are 19mm dia. and mounting 4 discs together is recommended for best performance.

Features & Benefits

- Three-dimensional abrasive/flexible bristle design
- Provide the convenience of a pre-charged abrasive but do not generate the heat associated with rubber wheels
- Coarse Disc for heavy cleaning
- Medium Disc for removing light oxide
- Fine Disk for pre-polishing, semi-finishing, blending, texturizing and polishing



keystone[®]
industries

dental.keystoneindustries.com

US: 480 South Democrat Road Gibbstown NJ, 08027
EU: Keystone Industries GmbH Stockholzstr. 11 78224 Singen, Germany

US: 1.800.333.3131
EU: +49 77 31 91 21 01

Applications

Orthodontic Appliances

Easy to smooth and polish orthodontic appliances without the need for emery paper.



Thermoplastic Resin

All thermoplastic resins, as well as nylon, polyester and polycarbonate, can be polished with HABRAS Disc Pro.



Metal Framework

Stippled plates can easily be cleaned.



COARSE YELLOW - 80 grit

Remove deep scratches from precious metals, non-precious metals, acrylic dentures and composites.

1670083, 4 Discs with Mandrel
1670097, 4 Discs without Mandrel
1670101, 48 Discs without Mandrel Bulk



MEDIUM RED - 220 grit

Removes fine scratches from precious metals, non-precious metals, acrylic dentures and composites.

1670084, 4 Discs with Mandrel
1670098, 4 Discs without Mandrel
1670102, 48 Discs without Mandrel Bulk



FINE BLUE - 400 grit

Smooths precious metals, non-precious metals, acrylic dentures and composites.

1670113, 4 Discs with Mandrel
1670085, 4 Discs without Mandrel
1670103, 48 Discs without Mandrel Bulk



SUPER FINE PINK - Pumice

Shine precious metals, non-precious metals, acrylic dentures and composites.

1670086, 4 Discs with Mandrel
1670100, 4 Discs without Mandrel
1670104, 48 Discs without Mandrel Bulk



HIGH GLOSS GREEN - 6μ

High Gloss precious metals, non-precious metals, acrylic dentures and composites.

1670108, 4 Discs with Mandrel
1670112, 4 Discs without Mandrel
1670116, 48 Discs without Mandrel Bulk



HIGH SHINE ORANGE - 1μ

High Shine precious metals, non-precious metals, acrylic dentures and composites.

1670107, 4 Discs with Mandrel
167011, 4 Discs without Mandrel
1670115, 48 Discs without Mandrel Bulk



HABRAS DISC KIT 1670113, 6 Discs with Mandrel

A unique three-dimensional, abrasive and flexible design that provides a uniform finish. Ideal for heavy cleaning, removing light oxide, pre-polishing, semi-finishing, blending, texturizing and polishing.

1 - Coarse Yellow - 80 grit 1 - Fine Blue - 400 grit 1 - High Gloss Green - 6μ
1 - Medium Red - 220 grit 1 - Super Fine Pink - Pumice 1 - High Shine Orange - 1μ

keystone[®]
industries

dental.keystoneindustries.com

US: 480 South Democrat Road Gibbstown NJ, 08027
EU: Keystone Industries GmbH Stockholzstr. 11 78224 Singen, Germany

US: 1.800.333.3131
EU: +49 77 31 91 21 01