

Trim[®] and Trim[®] II

Trim[®] and Trim[®] II are leaders in temporary crown and bridge. It is a fast-setting, non-irritating PEMA formulated for extended use.

Trim and Trim II are fast-setting, non-irritating PEMA temporary resin acrylic crown and bridge materials formulated for extended use. The powder and liquid in each kit are identical in formula and the powders can be mixed to obtain specialty colors. The resin is dimensionally stable, will not distort when withdrawn from undercuts while setting and is suitable for use as a bite registration material. Trim is available in two shades (light and dark). Trim II is available in six shades for greater diversity in achieving the perfect shade match for patients.

Recommended for use with:

- Glaze[®] Self Cure Resin Polish
- Masque[®] Silicone Separator & Lubricant
- Pressure Pot Hydraulic Water Press



Trim[®] Standard Kit (2-3 oz. powders, 6 oz. liquid)

PRODUCT	SIZE	ITEM NUMBER
Trim [®] Standard Kit		0921900
Light Powder	3 oz. (85g)	0921901
Dark Powder	3 oz. (85g)	0921902
Light Powder	1 lb. (454g)	0921907
Dark Powder	1 lb. (454g)	0921908
Liquid	4 oz. (118ml)	0921091
Liquid	6 oz. (177ml)	0921903
Liquid	8 oz. (236ml)	0921909
Liquid	32 oz. (946ml)	0921914



Trim[®] II Standard Kit (6 - 1.5oz. powders, 4 oz. liquid)

PRODUCT	SIZE	ITEM NUMBER
Powder, Shade B1 (59)	1.5 oz. (42g)	0921092
Powder, Shade A1 (60)	1.5 oz. (42g)	0921105-60
Powder, Shade A3 (66)	1.5 oz. (42g)	0921105-66
Powder, Shade B2 (61)	1.5 oz. (42g)	0921093
Powder, Shade A2 (62)	1.5 oz. (42g)	0921094
Powder, Shade D3 (65)	1.5 oz. (42g)	0921095
Powder, Shade C2 (69)	1.5 oz. (42g)	0921096
Powder, Shade B4 (77)	1.5 oz. (42g)	0921097
Bulk Powder, Shade B1 (59)	1 lb. (454g)	0921100-59
Bulk Powder, Shade A2 (62)	1 lb. (454g)	0921100-62
Bulk Powder, Shade D3 (65)	1 lb. (454g)	0921100-65
Bulk Powder, Shade C2 (69)	1 lb. (454g)	0921100-69
Bulk Powder, Shade B4 (77)	1 lb. (454g)	0921100-77
3ml Pipette Dropper	15	0921089



dental.keystoneindustries.com