

TruRepair®



TruRepair® is a Fast Bonding Acrylic Repair Material that is a must have for every office!

About TruRepair®

TruRepair is a fast-bonding acrylic repair material for repairing fractured dentures, dislodged denture teeth, and any additional office or lab emergency. A special internally-pigmented spherical globe of powder will blend with denture base materials to produce a color-stable match. TruRepair's self-curing cross-linked monomer and a PMMA powder system produces a high-strength repair and bond.

Features & Benefits

- Powder available in pink, vein and clear
- In many cases, a matrix will be needed. The use of the Bosworth Pressure Pot will increase the strength and the performance of the repair
- Mixing ratio: 2-3 parts powder to 1 part liquid



PRODUCT	SIZE	ITEM NUMBER
Liquid Only	4 oz.(236ml)	092164
Powder Only, Clear	3 oz. (85g)	092163
Powder Only, Vein	3 oz. (85g)	092162
Powder Only, Pink	3 oz. (85g)	092161
Standard Kit	 2 - 3 oz. (85g) Powder  4 oz. (118ml) Liquid	092160

keystone[®]
INDUSTRIES

dental.keystoneindustries.com

US: 480 South Democrat Road Gibbstown NJ, 08027
EU: Keystone Industries GmbH Stockholzstr. 11 78224 Singen, Germany

US: 1.800.333.3131
EU: +49 77 31 91 21 01