



## Patient is ready for a permanent reline when experience control and acceptance with TruSoft™

### About TruSoft™

TruSoft Resilient Denture Acrylic Relining Material is PEMA formulated and specifically designed as an intermediate liner for patients requiring a prolonged adjustment period prior to using a hard liner or denture base material. Its elasticity eliminates pressure pain, and accommodates tender ridges, sensitive areas, and irritated tuberosities. TruSoft is formulated to remain resilient for up to 30 days.

TruSoft can be used for any patient as a precautionary measure, but it is especially necessary for patients who are slow to heal or sensitive to wearing a new denture.



### Features & Benefits

- Bonding liquid strengthens adhesion to the denture
- Formulated to remain resilient for up to 30 days
- Elasticity ensures against pressure pain
- Gentle on tender ridges and sensitive areas

PRODUCT	SIZE	ITEM NUMBER
Standard Kit	Pink Powder	0921250
Powder Only	3 oz. (85 g)	0921251
Liquid Only	4 oz. (118 ml)	0921252
Bonding Liquid	1/3 oz. (10 ml)	0921983
Bulk Kit	1 lb. (454 g)	0921253
Powder Only	1 lb. (454 g)	0921254
Liquid Only	8 oz. (236 ml)	0921255
Bonding Liquid	4 oz. (118 ml)	0921984