

# Fastray™

## Custom Tray and Acrylic Base Plate Material

This classic product earned its' reputation as a leading custom tray material worldwide. Fastray™ is a dimensionally stable, radiopaque, self-curing individual tray and acrylic base plate material available in regular (5 min.) or extra-fast (3 min.) set. Fastray™ has no memory and when used with a separator it will not stick to fingers or the model. Available in 3 colors: Blue, Pink and White for different procedures.



### ORDERING INFORMATION

#### CUSTOM TRAY AND ACRYLIC BASE PLATE MATERIAL

- 0921375 Standard Kit, Regular Blue
- 0921376 Standard Kit, Regular Pink
- 0921377 Standard Kit, Regular White
- 0921378 Standard Kit, X-Fast Blue
- 0921379 Standard Kit, X-Fast Pink
- 0921380 Standard Kit, X-Fast White
- 0921381 Powder Only: Blue – 1 lb. (454g)
- 0921382 Powder Only: Pink – 1 lb. (454g)
- 0921383 Powder Only: White – 1 lb. (454g)
- 0921384 Regular Liquid Only – 8 oz. (236ml)
- 0921385 X-Fast Liquid Only – 8 oz. (236ml)
- 0921429 Powder \*Specify Color\* – 25lb (11.34kg)

- 0921386 Bulk Kit, Regular Blue
- 0921387 Bulk Kit, Regular Pink
- 0921388 Bulk Kit, Regular White
- 0921389 Bulk Kit, X-Fast Blue
- 0921390 Bulk Kit, X-Fast Pink
- 0921391 Bulk Kit, X-Fast White
- 0921392 Powder Only: Blue – 5lb (2.3kg)
- 0921393 Powder Only: Pink – 5lb (2.3kg)
- 0921394 Powder Only: White – 5lb (2.3kg)
- 0921395 Liquid Regular – 1 qt. (32oz., 946ml)
- 0921396 Liquid X-Fast – 1 qt. (32oz., 946ml)
- 0921415 Liquid Regular – 1 gal. (3.8l)
- 0921416 Liquid X-Fast – 1 gal. (3.8l)

# Fastray™

## LC Light-Cure Custom Tray Material

Fastray™ LC is a one component, light-cured custom tray material supplied in pre-formed wafer that, when used with a separator, will not distort or stick to fingers and/or models. Procedure time is reduced and odors are eliminated due to no mixing of powder and liquid.



### ORDERING INFORMATION

#### LC LIGHT-CURE CUSTOM TRAY MATERIAL

- 0921405 Standard Kit, Blue
- 0921406 Standard Kit, Pink
- 0921407 Standard Kit, White