

# Hatho Habras Miniature Brushes

Hatho Habras Miniature Brushes are manufactured in several grit sizes from long lasting polyimide filaments. All brushes are embedded with polish, eliminating the need for polishing paste.

The light grey brushes are embedded with aluminum oxide and are designed for pre-polishing ceramics and light cured materials. The dark grey brushes are embedded with silicone carbide and are designed for use on all alloys.

They last 12 times longer than standard goat hair brushes. The main advantages of Habras Brushes are durability, reduced dust while polishing, and a high sheen that can be achieved without the use of polishing pastes.



## ORDERING INFORMATION

### FINE

1670054 7/8" Aluminum Oxide  
1670055 1" Aluminum Oxide

### EXTRA FINE

1670068 3/4" Silicon Carbide  
1670069 7/8" Silicon Carbide  
1670071 1-5/8" Silicon Carbide  
1670072 2" Silicon Carbide

### MEDIUM

1670063 3/4" Aluminum Oxide  
1670064 7/8" Aluminum Oxide  
1670065 1" Aluminum Oxide  
1670073 3/4" Silicon Carbide  
1670074 7/8" Silicon Carbide  
1670076 1-5/8" Silicon Carbide  
1670077 2" Silicon Carbide

### COARSE

1670059 7/8" Aluminum Oxide  
1670060 1" Aluminum Oxide  
1670062 2" Aluminum Oxide Lathe

## UNIT SPECIFICATIONS

**SPEED** 3,000-6000 RPM  
**PACKAGE** 4 per package

## FEATURES AND BENEFITS

- 1 Durability & reduced dust while polishing
- 2 Lasts 12x longer than goat hair brushes
- 3 Brushes are embedded with polish
- 4 Several grit sizes ranging from: Fine, Extra Fine, Medium, & Coarse
- 5 Aluminum oxide made for pre-polishing
- 6 ceramics & light cured materials
- 7 Silicone carbide is made for alloys