


Original Truliner®

Original Truliner® is a color stable, permanent PMMA denture corrective relining material formulated for extended life (3-5 years).

ORDERING INFORMATION


ORIGINAL TRULINER®

- 0921940 Standard Kit, Pink
- 0921941 Standard Kit, Vein
- 0921942 Standard Kit, Clear



ORIGINAL TRULINER®

- 0921951 1 lb (454 g) Bulk Kit, Pink
- 0921952 1 lb (454 g) Bulk Kit, Vein
- 0921953 1 lb (454 g) Bulk Kit, Clear



FEATURES AND BENEFITS

- 1 Available in pink, vein and clear
- 2 Special monomer penetrates denture base, forming a chemical union & superior bond
- 3 Fine milled powder forms a dense relining preventing odors, tastes or stains
- 4 No-chip or peel formula

New Truliner®

New Truliner® is a PEMA denture corrective relining material designed to restore a denture to a patient's comfortable fit in a little as 10 minutes.

FEATURES AND BENEFITS


- 1 Available in pink, vein, and clear
- 2 Formulated to last up to 1 year
- 3 Color stable
- 4 No-sting, no burn formula



ORDERING INFORMATION

NEW TRULINER®

- | | | | |
|---------|---------------------|---------|------------------------------|
| 0921970 | Standard Kit, Pink | 0921977 | 1 lb (454 g) Bulk Kit, Pink |
| 0921971 | Standard Kit, Vein | 0921978 | 1 lb (454 g) Bulk Kit, Vein |
| 0921972 | Standard Kit, Clear | 0921979 | 1 lb (454 g) Bulk Kit, Clear |



TruRepair™

TruRepair™ is a fast-bonding, cross linked monomer fabricated to work with a fine milled acrylic repair material to repair fractured dentures, replace broken or dislodged denture

ORDERING INFORMATION

TRUREPAIR™

092160 Standard Kit (Specify Powder Shade)



2 - 3 oz. (85g)
Powder



4 oz. (118ml)
Liquid

RECOMMENDED FOR

USE WITH TRUREPAIR™

CORA-CAINE®

Denture Analgesic Adhesive

PRESSURE POT

Hydraulic Water Press

MASQUE®

Silicone Lubricant & Separator



FEATURES AND BENEFITS

- 1 Use for repairing fractured dentures or temporary partials
- 2 Available in pink, vein and clear
- 3 Color stable
- 4 Cross linked