

Tecnoflask



Tecnoflask allows Diamond D® and other acrylics to be processed in much quicker time inside of a microwave.

About Tecnoflask

The TF-100 muffle will give your laboratory a cleaner environment and make the process more enjoyable for the production of your work. Included with the kit and available separately is the PL-100 protection plate, which extends the muffle's service life with extra needed protection.

Features & Benefits

- Acrylic material designed for use in microwave
- Cannot be damaged by exertion of pressure
- Polymerization time shortened to four minutes
- Used in the making of total upper & lower dental prostheses
- No isolating material needed between the muffle & the plaster
- Does not require a specific resin

Recommended for use with:

» Diamond D® Denture Acrylic » Diamond D® Ultra Sep



PRODUCT	ITEM NUMBER
Tecnoflask Kit	1009198
Tecnoflask Muffle TF-100	1009200
Tecnoflask Press Board PL-100	1009199

keystone[®]
industries

dental.keystoneindustries.com

US: 480 South Democrat Road Gibbstown NJ, 08027
EU: Keystone Industries GmbH Stockholzstr. 11 78224 Singen, Germany

US: 1.800.333.3131
EU: +49 77 31 91 21 01