

# Machine III

## The Powerful Vacuum Former



The Machine III has been a staple in Keystone Industries' long-lasting line of highly recommended vacuum formers.

The price, features and reliability make it the first choice in dental vacuum formers. The machine has more power than any other dental vacuum former, which means the production of more accurate and definite dental thermoplastic appliance.

The aluminum construction is long-lasting, eliminating the possibility of rust. The heating element is symmetrically engineered for even heat distribution throughout the thermoplastic material.

This machine is used to process bleaching trays, sports mouthguards, implant stints, temporary splints, custom impression trays, denture bases, copings, temporary crowns and bridges and orthodontic appliances.

### RECOMMENDED FOR USE WITH



Pro-Form® Mouthguard Laminates



### FEATURES AND BENEFITS

- Easy-to-use system
- Uses 5"x 5" resin sheets
- Symmetrical heating system
- Lightweight, rust-free aluminum frame
- Made in the U.S.A.
- 5-year guarantee on parts and labor

### ORDERING INFORMATION

Machine III, 110 Volt, UL/cUL/CE marked	#7000330
Machine III, 220 Volt, UL/cUL/CE marked	#7000332

