

Pro-Form[®] Mouthguard Laminates

Dentistry's Leader in Custom Laminates

Pro-Form manufactures the only pre-made laminated mouthguard material in the industry.

This material has long been considered the leader in custom-made mouthguards. By laminating two layers of soft EVA together, the laminate has excellent tensile strength. All of these features make the Pro-Form mouthguard laminate the material of choice for custom-made mouthguards.



FEATURES AND BENEFITS

- Available in .160 (4mm) 5"x5" (127 x 127mm) squares
- Also available in round
- Available in 12/Pkg
- Available in assorted color packs

SINGLE COLOR



-  #9601020 White
-  #9601040 Yellow
-  #9600850 Orange
-  #9600910 Red
-  #9600886 Pink
-  #9600820 Maroon
-  #9600870 Purple
-  #9598060 Blue
-  #9600800 Green
-  #9598730 Black
-  #9600850 Clear

DUAL-COLORED



TRI-COLORED



TIE-DIE GUARDS



GLITTER GUARDS



FUN GUARDS



-  #7924000 Wilderness
-  #7954000 Patriot
-  #7944000 Midnight
-  #7965000 Shamrock
-  #7934000 Flamingo
-  #7911111 Rainbow
- #7974000 All In

-  #7965030 Yellow Glitter
-  #7965020 Orange Glitter
-  #7965050 Purple Glitter
-  #7965040 Green Glitter
-  #7965060 Pink Glitter
-  #7965010 Blue Glitter
- #7965070 Assorted

-  #9601039 White Polkadot
-  #9599077 Blood
-  #9598079 Blue Polkadot
-  #9599075 Fangs
-  #9600929 Red Polkadot

