

# ProphyBrite® Air Polisher

## Brite Smiles All Around

With this dental air polisher, clogging and wrist fatigue will be the least of any dental hygienist's worries.

It has a unique system and ergonomical design. The ProphyBrite provides quick and gentle patient care and requires no installation. It works best by delivering a fine sodium bicarbonate and water slurry for easy removal of plaque and stains. The powder flows through this handpiece smoothly and efficiently, ensuring it cleans only the desired areas of the mouth. There are very minimal clogs when using this unit, as long as it gets cleaned properly. It's the ideal dental air polisher for cleaning and preparing teeth prior to placement of orthodontic brackets, sealants, bleaching, and fluoride treatments.



### FEATURES AND BENEFITS

- Lightweight & Ergonomic
- Autoclavable Nozzle
- Large Powder Chamber (13G)
- Available in orange, black, red, blue
- No installation needed
- Rotating Nozzle
- Large Variety of Adapters

### RECOMMENDED FOR USE WITH



Prophy Powder

### KITS

|                                 |          |
|---------------------------------|----------|
| Air Polishing Unit (Red)        | #16670R  |
| Air Polishing Unit (Black)      | #16670B  |
| Air Polishing Unit Kit (Orange) | #16670N  |
| Air Polishing Unit Kit (Blue)   | #16670LB |

### ADAPTORS

|                       |        |
|-----------------------|--------|
| Adaptor - Kavo        | #16678 |
| Adaptor - w&H-Adec    | #16684 |
| Adaptor - NSK Mach-QD | #16680 |
| Adaptor - NSK Ti-Max  | #16683 |
| Adaptor - Bien Air    | #16671 |
| Adaptor - Dental EZ   | #16686 |

### COMPONENTS

|   |         |
|---|---------|
| Nozzle  | #16677  |
| Cleaning File Set (1-2.5cm, 1-8cm, 1-10cm)                | #16674  |
| Spare Parts Set (1 powder Chamber Cap & 1 set of O-rings) | #16681K |