

Versacryl®

Not Your Average Soft Reline

This unique product has been developed to create thermo-adjustable relines.

Imagine having the luxury to simply heat a denture and re-adjust the fit in seconds and if it's not satisfactory the procedure can be repeated anytime by the dentist, the technician or even the patient themselves.

The Versacryl Reline is not a regular "soft reline" material. It is a heat sensitive liner that softens in warm water between 120F – 205F. At body temperature the Versacryl has an extremely comfortable cushiony effect that yields to pressure and actually creates better suction than conventional soft liners.

Versacryl is a cross-linked acrylic denture material guaranteed to create a chemical bond with all other denture acrylics. Versacryl is available in Heat Cure or Self Cure and no special equipment is needed to process.



Versacryl Internal Reline
Gasket partial/Itsoclear clasps

Versacryl Gasket Retention Partial

FEATURES AND BENEFITS

- Full and Partial denture relines
- Added Lingual retention on lower dentures
- Obturator cases (cancer patients)
- Post dam areas on full uppers
- Gasket retention partial dentures
- Lingual retention on partial dentures
- Relining flexible partials
- Pedio Partial cases
- Retention under acrylic clasps
- Internal clear soft layers on nightguards and bite splints

ORDERING INFORMATION

Versacryl Heat Cure Lab Kit	#1014032
Versacryl Heat Cure Hardener Liquid – 60ml	#1014012
Versacryl Heat Cure Pre-Mixed – 60ml	#1014014
Versacryl Heat Cure Hardener – 1 qt.	#1014017
Versacryl Heat Cure Pre-Mixed – 1 qt.	#1014019
Versacryl Heat Cure Reline Kit	#1014034
Versacryl Softening Liquid – 60ml	#1014010
Versacryl Softening Liquid – 1 qt.	#1014015
Versacryl Bulk Clear Powder – 1 lb.	#1014022
Versacryl Bulk Pink Powder – 1 lb.	#1014023
Versacryl Clear Powder – 2 oz.	#1014020
Versacryl Pink Powder – 2 oz.	#1014021
Versacryl Sealer – 1/2 oz.	#1014024