

Support Peg

Instant firing support

Support Peg is one of the first of its kind, designed specifically to simplify the dental firing process.

It offers instant firing support by holding the crown or jacket in place.

To use, simply fill the crown or jacket and insert the Support Peg syringe; no drying time is necessary. Support Peg will hold the crown in place for the entire firing process, and once complete, will crumble out with ease proving to be time saving yet efficient. Unlike other methods, Support Peg fires without any shrinkage or distortion of the structure when inside the furnace.



FEATURES AND BENEFITS

- For use with e.max and Lava crowns
- No pre-drying time
- Zero firing distortion
- Does not contaminate
- For temperatures from 400°C – 1300°C

ORDERING INFORMATION

Support Peg

#1009301

