

# Clearmet

## The Clear Alternative

NEW

Clearmet is a revolutionary new clear thermoplastic resin specifically designed for durability and biocompatibility.

The material is intended for use in an injection molding process. Unlike conventional metal partials, Clearmet is an ultra-transparent framework. Patients can now enjoy a lightweight partial without a metal taste or undesirable metal clasps showing on their teeth. No one will ever know a partial denture restoration is being worn.

Clearmet is made of monomer-free material that has no known allergic reactions. The nonporous, bacteria and stain resistant material means there will be no odors given off as well.

The goal of Clearmet is not only to improve smiles, but more importantly, to improve lives. Patients will be more confident than ever to flash their new appearance! Clearmet is the only “clear choice” for denture partials.



### FEATURES AND BENEFITS

- Frame is invisible in the mouth
- Stain resistant, odor-free material
- Can be relined & repaired in saddle areas
- Clasps are tooth bearing
- 90° into undercut creates excellent retention
- Clasps are adjustable with a warm three-pong plier, making it easy to adjust chairside
- Can be used with Myerson, TCS, & Valplast machines

### ORDERING INFORMATION (5 tubes per package)

Clearmet, Small Tubes (1.77 in.)	#1009225
Clearmet, Medium Tubes (3.03 in.)	#1009227

### RECOMMENDED FOR USE WITH



Diamond D® Denture Acrylic



Diamond D® Ultra Sep

### IMPORTANT NOTE:

Prior to injection process, follow instructions on necessary moisture removal.



The Clearmet denture partial doesn't leave any unsightly metal clasps, unlike other denture partials. The frame is invisible in the mouth.

keystone<sup>®</sup>  
I N D U S T R I E S



keystone  
industries

480 South Democrat Rd. Gibbstown, NJ

## About Us: The Keystone Mission

Since 1908, Keystone Industries has maintained a strong reputation of producing innovative high-tech dental products in both the operator and laboratory realms. We are committed to the goal of total customer satisfaction. We are dedicated to the achievement of this goal by meeting or exceeding our customer expectations for product quality, competitive pricing and prompt delivery and customer service.



### Video Learning Center

Our new "In The Know" video series gives you an up close look at hands-on demonstrations of products being done by professionals in the field.



### Product Literature

Download all catalogs, brochures, booklets, sell sheets, SDS, and instructions from our online databases.




### Promotions


We always have our quarterly promotions posted on our website. Check them out and take advantage of the best dental deals in the industry.


## Contact Information

### Domestic Office

480 South Democrat Road  
Gibbstown NJ 08027

 800-333-3131


 856-663-4700


 856-224-9444

 [customerservice@keystoneind.com](mailto:customerservice@keystoneind.com)

### International Office

Werner-von-Siemens Str. 14a,  
D-78224 Singen, Germany

 +49 77 31 91 21 01

 +49 77 31 91 21 02

 [sales@keystoneind.com](mailto:sales@keystoneind.com)