


BEST
SELLER

True refinement
requires

Hatho Miniature Scotch Brite™



**BUY any 3 dozen
Hatho Miniature Scotch Brites™,
GET 1 at no charge!**

Hatho Miniature Scotch Brite™

You can use the Hatho Miniature Scotch Brite™ wheel to polish services without ruining delicate anatomical surfaces on dental appliances. By using low speeds – around 5,000 rpm – and light pressure the acrylic should not overheat due to friction. These brushes come in coarse, medium and fine to let you control the type of polishing finish you need.

- Available in 12/Pkg
- Available in Coarse, Medium, Fine



Coarse
(#1670016)



Medium
(#1670092)



Fine
(#1670091)

This offer cannot be combined with any other promotion or discount. Other restrictions may apply. Customer will be invoiced for purchased goods from preferred distributor and no charge goods will ship with the same order. Offer valid from 01/01/2019—04/09/2019. Valid in the US only. Not responsible for typographical errors.

TO ORDER, CALL: 1.800.333.3131 OR VISIT DENTAL.KEYSTONEINDUSTRIES.COM