

ProCure 300 Light Oven

The ProCure 300 Light Oven is an advanced third generation light curing unit that provides superior lighting intensity and excellent depth of cure, without the negative heat buildup.

- A highly reflective chamber for uniform light coverage
- A rotating glass platform to ensure uniform light
- A built-in fan to maintain optimum temperature to eliminate warping
- Two powerful halogen bulbs to ensure excellent intensity and depth of cure
- A protective welder-quality shield of glass



**BUY a ProCure Oven,
GET a Bosworth Light Liner[®]
at no charge!**

#9634590

ProCure 300 Light Oven

#9634600

ProCure 300 Light Oven with sample pack



**BUY a 12/pk Pro-Form Laminates,
GET a Scotch Brite™ at no charge!**

See website for more information

Pro-Form[®] Mouthguard Laminates

Our Pro-Form Mouthguard Laminates are the only pre-made laminated mouthguard material in the industry, and are widely used by leaders in custom-made mouthguards. By laminating two layers of soft EVA together, the laminate has excellent tensile strength. All of these features make the Pro-Form mouthguard laminate the material of choice for custom-made mouthguards.

- Available in .160 (4mm) 5"x5" (127 x 127mm) squares, as well as 120mm or 125mm round laminates
- 12/pack
- Colors Available: Clear, Orange, Yellow, Maroon, White, Green, Black, Red, Blue, Purple, Pink, as well as Assorted color packs

Pro-Form[®] Electric Knife

The Pro-Form Electric Knife is useful in trimming and cutting a variety of thermoplastic dental materials. The unit has solid state circuitry, an interchangeable handpiece, 2 outlet stations to accommodate more than one person, all stainless steel construction, and 10 amps of power.

- 10 amps allow you to cut most materials easily and smoothly by providing extra heat at the tip
- Excellent for bleaching trays, impression trays, splints, mouthguards and more.



**BUY an Electric Knife,
GET 2 12/pk of Pro-Form Laminates
at no charge!**

#9144730

Pro-Form Electric Knife

This offer cannot be combined with any other promotion or discount. Other restrictions may apply. Customer will be invoiced for purchased goods from preferred distributor and no charge goods will ship with the same order. Offer valid from 04/01/2019—06/30/2019. Valid in the US only. Not responsible for typographical errors.

FOR MORE INFORMATION, VISIT DENTAL.KEYSTONEINDUSTRIES.COM. TO ORDER, CALL YOUR PREFERRED DISTRIBUTOR.