

+OPTIMIZE

ORANGE SOLVENT

ORANGE SOLVENT WITH LANOLIN


Used in the removal of zinc oxide-eugenol,
rubber base and other impression materials
from hands, face, mixing slabs and instruments.



Keep out of reach of children.

Net Contents:
32 fl oz (946 ml)

REF 00134



Directions:

Moisten cloth or swab with **Orange Solvent** and wipe away
impression material.

Orange Solvent cleans quickly, easily and is non-irritating.



Caution: Do not take internally.

Contains lanolin. Those sensitive to lanolin should not use this product.

Warning

Causes serious eye irritation. May cause an allergic skin reaction. Refer
to product SDS for additional health, safety and disposal information.



Keystone Industries
52 West King St.
Myerstown, PA 17067 USA
800.333.3131
keystoneindustries.com

Made in USA

keystoneindustries.com/symbolglossary

International Sales

Keystone Industries GmbH,
Stockholzstr. 11, 78224 Singen, Germany
+49 77 31 91 21 01

EPR-18 NatLab-1364 4926630 02/2021