

PRECISION 3D RESINS

KeySplint Soft™



#4220005 KeySplint Soft, 35.3 oz (1 kg)

FEATURES OF KEYSPLINT SOFT

KeySplint Soft is fully biocompatible, and is in compliance with international medical device regulations and standards. The robust formulation has two years of guaranteed shelf life, is color stable, and is validated in multiple open-source printers. No other dental splint resin can deliver the quality and performance of KeySplint Soft.

510k-Approved, EU Class IIa Approved

- Biocompatible
- Transparent
- Polishable
- Stain-resistant

PROPERTY	PROCEDURE	KEYSPLINT SOFT
Flexural Strength	ASTM D790	47
Flexural Modulus	ASTM D790	1356
Tensile Strength	ASTM D638	52
Young's Modulus Tensile	ASTM D638	1790
Water Sorption	ISO 20795-2	<84.8 µg/mm ³
Biocompatibility	ISO 10993-1 and USP Class VI	Pass