



Tray Adhesive Compound (TAC)

TAC™ Tray Adhesive Compound is a leading tray adhesive compound available in convenient spray cans or brush-on liquid for use with alginate impression materials.

About TAC™

TAC Tray Adhesive Compound is a leading tray adhesive compound available in convenient spray cans or brush-on liquid for use with alginate impression materials. TAC Tray Adhesive Compound securely holds the alginate in the tray when the impression is removed from the patient's mouth. TAC Tray Adhesive Compound is alcohol-based for faster drying and easier clean-up and can be dissolved with Isopropyl Alcohol.

Features & Benefits

- Securely holds the alginate in the tray when the impression is removed from the patient's mouth
- Alcohol-based for faster drying
- Easy clean-up
- Can be dissolved with Isopropyl Alcohol
- Can be used as an adhesive to hold Softone™ and Trusoft™



PRODUCT	SIZE	ITEM NUMBER
Standard Kit	2 x .5 oz (14.8 ml) Bottles (Brush-on Liquid)	0921881
Bulk Bottle	4oz (118ml)	0921882
Bulk Bottle	8oz (236ml)	0921883
Spray	3.5oz (99g)	0921880



dental.keystoneindustries.com