



## SUPER-VAC 40 PLUS

Latest generation of concentrated evacuation cleaner. Powerful multi-enzyme cleaner that deodorizes as it cleans with a pleasant mint fragrance. EPA-Compliant and safe for use with amalgam separators. The non-foaming, environmentally friendly, neutral pH formula dissolves organic fluids and debris. Recommended for daily use to maintain efficacy of evacuation equipment. Makes 8 gallons of solution.



PRODUCT	SIZE	ITEM NUMBER
Liquid Concentrate	64 fl oz (1.89 L)	00138

Pairs well with our HVEs and Saliva Ejectors.



[dental.keystoneindustries.com](http://dental.keystoneindustries.com)