



TEMPORARY CEMENT REMOVER

This solution is ready to use as is. This strong fast-acting formula removes cement and stone quickly and completely from instruments, trays and other appliances. Superior cleaner for most zinc-oxide eugenol type cement.

PRODUCT	SIZE	ITEM NUMBER
Liquid	128 fl oz (3.78 L)	00140



dental.keystoneindustries.com

US: 480 South Democrat Road Gibbstown NJ, 08027
EU: Keystone Industries GmbH Stockholzstr. 11 78224 Singen, Germany

US: 1.800.333.3131
EU: +49 77 31 91 21 01